


Project Pipeline



\$ Rates

As of July 2018 Users no longer need to bring enzymes for their runs. Enzyme rates have been incorporated into our rates.


Project set up fee: \$782

Example: 4 samples in duplicate \$3223

See SCC website for full rate schedule

* sequencing costs not covered in SCC fees

Mandovi Chatterjee, Alex Ratner, Ian Brenckle, Sarah Boswell



InDrops method overview: cells are encapsulated into droplets with lysis buffer, reverse-transcription mix, and hydrogel microspheres carrying barcoded primers. After encapsulation primers are released. cDNA in each droplet is tagged with a barcode during reverse transcription. Droplets are then broken and material from all cells are linearly amplified before sequencing.


How to Succeed with scRNA-Seq?

- 1) Practice sample dissociation prior to run date.
 - a) Gentle cell dissociation: shorter time, lower temperature.
 - b) OptiPrep density gradient: cushion cells to avoid pelleting.
 - c) Try MACS sorting instead of FACS for increased viability.
- 2) Ensure cells numbers and viability >90% viable after 1 hour on ice.
 - a) Strain cells & confirm cells do not start clumping.
 - b) Count cells on hemocytometer. Need 15,000 cells minimum.
- 3) Work with SCC to optimize experimental design to address:
 - a) Batch effects & number of replicates
 - b) Number of single cells to profile
- 4) Bring cells in buffer WITHOUT calcium, EDTA, or heparin as these will inhibit RT. Magnesium levels should be <3mM.
 - a) Suggest PBS with 0.1-1% BSA

Sample Handling

- Cells kept on ice until lysis in droplet
- RT occurs within single cell droplet
- Samples frozen until library preparation

Microfluidic Device




Library Preparation

- Clean up of unused primers from DNA:RNA hybrid
- Second strand synthesis to full length cDNA
- In vitro transcription to amplify product as RNA
- Fragmentation and RT by random hexamers
- PCR adaptor addition and library amplification

Data Analysis

- Experienced informaticists can follow analysis code.
- 10x users can use Cell Ranger Software.
- Suggest work with Chan Bioinformatics Core <http://bioinformatics.sph.harvard.edu/>

New 10x Chromium Single Cell Controller New



10x Chromium scRNA-Seq is available as a staff assisted or self-service model. Talk to the SCC to learn more!